



Better Health. Better Future.

Phlebotomy Blood Drawing Chairs

Quality and easy-to-use phlebotomy solutions to improve patient care.



Experience the difference of a Tech-Med® blood draw chair

Ease-of-Use

The cylindrical design and one-control knob allow users to easily raise, lower, and rotate the full length arms into infinite positions, ensuring optimal setup and comfort for both phlebotomist and patient.

Patient Comfort

Adjustable flip arm can extend out up to 7.5" to comfortably accommodate patients of all sizes, as well as, a child sitting on a parent's lap. 7" wide arm allows for patient comfort and greater flexibility for the phlebotomist.

Durability

The welded solid steel frame, 300-700 lb weight capacity, and adjustable levelers on all four legs provide long-lasting durability, strength, and stability.



4383-F
Upholstered, Standard Height
Weight capacity: 500 lb



4382F-BLK
Upholstered Standard Height
Weight capacity: 500 lb

Upholstered blood draw chairs feature tool-free assembly, making setup as easy as 1-2-3



Easy-to-use, one-control knob accommodates patient and phlebotomist needs with infinite arm positions



Available in a variety of colors



Adjustable floor levelers allow chair to sit solidly on uneven surfaces



4382-F
Upholstered Standard Height
Weight capacity: 500 lb



One-piece adjustable flip arm



Expands up to 7.5" for a patient's comfort

Optimize Your Space With Our Foldable Design

Tech-Med's® foldable blood draw chair is a perfect fit for mobile, temporary, or space-limited areas. Engineered with comfort and convenience in mind, it features plush padding for patient comfort and a foldable design for easy transport and storage. Its solid iron construction and maneuverable frame ensure swift setup and repositioning, instilling confidence in efficient blood collection in any setting.



Chair folded

Product dimensions closed:
26" width, 6" depth, 44" height

No assembly required

4380
Blood Draw Chair, Folding
Product dimensions:
26"width x 22"depth x 37"height
Weight capacity 500lb



Tech-Med® tray and drawer accessory

Keep phlebotomy supplies on hand and within reach with our tray and drawer accessory. The pre-drilled straight arm of Tech-Med blood draw chairs ensures easy installation on all models.

4388
Tray & Drawer for Blood Draw Chair

Tech-Med® Lab Blood Draw Chair

Flip arm expands up to 7.5" for patients of all sizes



One control knob to easily raise and lower



Rotate the full length arm into infinite positions



Adjustable floor levelers allow chair to sit solidly on uneven surfaces

4385F-BLK
Lab Series Tall
Weight capacity: 300 lb



All parts, tools, and instructions included

The four-bolt design of Tech-Med lab chairs ensures stability and strength in an easy to assemble design. Each of the four bolts secure the seat to the frame allowing one individual to complete assembly easily and quickly.



Leading National Brand of Products

Our extensive line of medical supplies provides quality and easy-to-use solutions to serve the everyday demands of healthcare professionals, with the belief that enhancing family and community health is the way to a better future.

Combined with the highest standards in supply chain management and service. We strive to be a trusted partner, so healthcare professionals can focus solely on the well-being of their patients.



Exam room supplies

Medical equipment that works as hard as you. Dukal provides a wide range of durable exam room supplies to provide your patients with the highest level of care.



Patient disposables

Your partner in delivering care, Dukal's comprehensive portfolio provides quality disposable solutions, including tourniquets, bandages, cleansing and antiseptic pads, and a wide range of white goods to meet the daily needs of phlebotomists.



Alcohol Pads



Sundry Jars



Design Series Adhesive Bandages



Cloth Tape



Tourniquets



Disposables

Tech-Med® Blood Drawing Chairs

Upholstered

Provide ultimate comfort to patients with padded seat, back, and arm.

Standard Tall Weight capacity **500 lb** ✓ Tool free assembly

Item #	Description	Overall Dimensions
4382	Standard Height, Straight Arm, Blue	34" width, 28" depth, 36" height
4382-F	Standard Height, Flip Arm, Blue	34" width, 28" depth, 36" height
4382F-BLK	Standard Height, Flip Arm, Black	34" width, 28" depth, 36" height
4383-F	Tall Height, Flip Arm, Blue	34" width, 28" depth, 43" height
4383F-BLK	Tall Height, Flip Arm, Black	34" width, 28" depth, 43" height



Upholstered Wide Bariatric

Extra room for individual patients and added comfort for younger patients with enough room for adult and child.

Standard Tall Weight capacity **700 lb** ✓ Tool free assembly

Item #	Description	Overall Dimensions
4382X-F	Standard Height, Flip Arm, Blue	43" width, 28" depth, 40" height
4382XF-BLK	Standard Height, Flip Arm, Black	43" width, 28" depth, 40" height
4386X-F	Tall Height, Flip Arm, Blue	43" width, 28" depth, 43" height
4386XF-BLK	Tall Height, Flip Arm, Black	43" width, 28" depth, 43" height



Lab Series

Available with upholstered arms and plastic seat to meet the facilities needs.

Standard Tall Weight capacity **300 lb**

Item #	Description	Overall Dimensions
4384F-BLK	Standard Height, Upholstered Flip Arm, Black	34" width, 28" depth, 32" height
4385F-BLK	Tall Height, Upholstered Flip Arm, Black	34" width, 28" depth, 43" height



Foldable

Folding design offers ease of transport and storage. No assembly.

Standard Weight capacity **500 lb**

Item #	Description	Overall Dimensions
4380	Blood Draw Chair, Folding	26" width, 22" depth, 37" height



2-year warranty on all Dukal blood draw chairs

> Find our complete product specifications at dukal.com



Dukal, LLC
2 Fleetwood Ct, Ronkonkoma, NY 11779

Phone: 1 631 656 3800
Toll Free: 1 800 243 0741

Email: customerservice@dukal.com
To Place Orders: orders@dukal.com

Find our full line of products at [dukal.com](https://www.dukal.com)