



## Product Update: Gauze Sponges

We are transitioning our packaging to reflect Dukal's recent rebrand. Below are the products that are affected. You will begin to see this update in the coming months.

### Non-Sterile Rolled Gauze

- 2-Ply 100% woven cotton gauze
- Absorbent, effectively holds dressings in place



Item #	Description	Packaging
403	3" x 5 yd 2-Ply	12/BG, 8 BG/CS

Item #	Description	Packaging
404	4" x 5 yd 2-Ply	12/BG, 8 BG/CS

### Sterile Rolled Gauze

- 2-Ply 100% woven cotton gauze
- Absorbent, effectively holds dressings in place



Item #	Description	Packaging
453	3" x 5 yd 2-Ply	1/PK, 12 PK/BX, 8 BX/CS

Item #	Description	Packaging
454	4" x 5 yd 2-Ply	1/PK, 12 PK/BX, 8 BX/CS

Find our full line of products at [dukal.com](http://dukal.com)

Dukal, LLC  
2 Fleetwood Ct, Ronkonkoma, NY 11779

Phone: 1 631 656 3800  
Toll Free: 1 800 243 0741

Email: [customerservice@dukal.com](mailto:customerservice@dukal.com)  
To Place Orders: [orders@dukal.com](mailto:orders@dukal.com)



## Product Update: Gauze Sponges

We are transitioning our packaging to reflect Dukal's recent rebrand. Below are the products that are affected. You will begin to see this update in the coming months.

### Sterile X-Ray Detectable Type VII Gauze Sponges

- USP Type VII gauze, 100% woven-cotton
- Blue filament runs through each sponge and has a unique wave pattern to ensure easy x-ray detection
- Sterile sponges are packaged in rigid plastic tray with peel back lid
- Both sterile and non-sterile sponges are banded in 10s to ensure accurate count



Item #	Description	Packaging
416-2x	4" x 4" 16-Ply	10/TY, 128 TY/CS

### Non-Sterile Type VII Gauze Sponges

- USP Type VII gauze, 100% woven cotton
- Highly absorbent and offers exceptional clinical performance
- Used for a number of applications including: general wound care, dressings, prepping, wound packing, and debriding



Item #	Description	Packaging
C83082	3" x 3" 12-Ply	200/BG, 20 BG/CS

Find our full line of products at [dukal.com](http://dukal.com)

Dukal, LLC  
2 Fleetwood Ct, Ronkonkoma, NY 11779

Phone: 1 631 656 3800  
Toll Free: 1 800 243 0741

Email: [customerservice@dukal.com](mailto:customerservice@dukal.com)  
To Place Orders: [orders@dukal.com](mailto:orders@dukal.com)