



Product Update: Sterile Gauze Sponges

We are transitioning our packaging to reflect Dukal's recent rebrand. Below are the products that are affected. You will begin to see this update in the coming months.

Sterile Gauze Sponges

- 100% woven cotton
- Used for a number of applications including: general wound care, dressing, prepping wound packing, and debriding
- Non-sterile available in 200 or 100 per bag to suit multiple market needs



Item #	Description	Packaging
6208	2" x 2", 8-ply, ST	2/PK, 50 PK/BX, 30 BX/CS



Item #	Description	Packaging
6312	3" x 3", 12-ply, ST	2/PK, 50 PK/BX, 24 BX/CS



Item #	Description	Packaging
6408	4" x 4", 8-ply, ST	2/PK, 50 PK/BX, 12 BX/CS



Item #	Description	Packaging
6842	8" x 4", 12-ply, ST	2/PK, 25 PK/BX, 20 BX/CS

Find our full line of products at dukal.com

Dukal, LLC
2 Fleetwood Ct, Ronkonkoma, NY 11779

Phone: 1 631 656 3800
Toll Free: 1 800 243 0741

Email: customerservice@dukal.com
To Place Orders: orders@dukal.com



Product Update: Sterile Gauze Sponges

We are transitioning our packaging to reflect Dukal's recent rebrand. Below are the products that are affected. You will begin to see this update in the coming months.

Product category update: Basic Care is now "Everyday Care"

Everyday Care Gauze Sponges

- 100% woven cotton
- Lighter, economical choice for quality care



Item #	Description	Packaging
8500	2" x 2", 8-ply, ST	2/PK, 50 PK/BX, 30 BX/CS



Item #	Description	Packaging
8502	4" x 4", 8-ply, ST	2/PK, 50 PK/BX, 12 BX/CS



Item #	Description	Packaging
8503	4" x 4", 12-ply, ST	2/PK, 25 PK/BX, 24 BX/CS



Item #	Description	Packaging
8505	2 x 2 8-Ply, NS	200/BG, 25 BG/CS

Find our full line of products at dukal.com

Dukal, LLC
2 Fleetwood Ct, Ronkonkoma, NY 11779

Phone: 1 631 656 3800
Toll Free: 1 800 243 0741

Email: customerservice@dukal.com
To Place Orders: orders@dukal.com



Product Update: Sterile Gauze Sponges

We are transitioning our packaging to reflect Dukal’s recent rebrand. Below are the products that are affected. You will begin to see this update in the coming months.

Product category update: Caliber™ is now “Premium Care”

Premium Care Type VII Gauze Sponges

- Meets USP Type VII requirements
- 100% woven cotton
- Highly absorbent and offers exceptional clinical performance
- Used for a number of applications including: general wound care, dressing, prepping wound packing, and debriding



Item #	Description	Packaging
C5119	2" x 2", 8-ply, ST	2/PK, 50 PK/BX, 30 BX/CS



Item #	Description	Packaging
C5139	4" x 4", 12-ply, ST	2/PK, 25 PK/BX, 24 BX/CS

Find our full line of products at dukal.com

Dukal, LLC
2 Fleetwood Ct, Ronkonkoma, NY 11779

Phone: 1 631 656 3800
Toll Free: 1 800 243 0741

Email: customerservice@dukal.com
To Place Orders: orders@dukal.com