



PASSION. PARTNERSHIP. POSSIBILITIES.



Caliber™ Surgical Drapes

Description

Tissue-poly-tissue drapes for general use in surgical procedures

Features & Benefits

- Polyethylene middle layer helps eliminate risk of bacterial migration and moisture strike through
- Use blue tissue side up to prevent glare
- Available in plain for general use and fenestrated to allow access to operative site while protecting the surrounding area
- Not made with natural rubber latex

Common Applications:

- Minor surgical procedures

Item#	Description	Packaging
Sterile		
20-001	Fenestrated, 18" x 26", 3-ply	1/pk, 50/bx, 6 bx/cs
20-002	Non-Fenestrated, 18" x 26", 3-ply	1/pk, 50/bx, 6 bx/cs

Connect with us!

Find our full line of products online at:

www.dukal.com

Telephone: 1-631-656-3800

Toll Free: 1-800-243-0741

Fax: 1-631-656-3810

Email: customerservice@dukal.com

To Place Orders: orders@dukal.com

Corporate Office

DUKAL Corporation
2 Fleetwood Court,
Ronkonkoma, NY 11779



Control # SL070719