



PASSION. PARTNERSHIP. POSSIBILITIES.



# Cotton Undercast Padding

## Description

Padding provides cushioning under casting material and skin; while absorbing moisture and enabling the skin to breathe

## Features & Benefits

- Constructed of 100% soft, cotton fibers
- Holds padding in place without shifting or bunching
- Conforms to all body contours
- Easy to tear

## Common Applications:

- Used with either synthetic or plaster casting material
- Padding splints

Item	Description	Packaging
Non-Sterile - Latex Free		
30122	2" x 12'	24 rl/bg, 3 bg/cs
30123	3" x 12'	12 rl/bg, 6 bg/cs
30124	4" x 12'	12 rl/bg, 6 bg/cs
30126	6" x 12'	6 rl/bg, 6 bg/cs

### Connect with us!

Find our full line of products online at:  
[www.dukal.com](http://www.dukal.com)

Telephone: 1-631-656-3800  
Toll Free: 1-800-243-0741  
Fax: 1-631-656-3810  
Email: [customerservice@dukal.com](mailto:customerservice@dukal.com)  
To Place Orders: [orders@dukal.com](mailto:orders@dukal.com)

### Corporate Office

DUKAL Corporation  
2 Fleetwood Court,  
Ronkonkoma, NY 11779

