



PASSION. PARTNERSHIP. POSSIBILITIES.



Xeroform Petrolatum Gauze

Description

Provides soothing non-adherent protection with more convenience, less mess, and freshens as it protects

Features & Benefits

- Non-adherent petrolatum
- Clings and conforms to all body contours
- Non-toxic, non-irritating, guaranteed sterile
- Contains 3% bismuth tribomophenate, a special deodorizing agent
- Not made with natural rubber latex

Common Applications:

- An initial layer in surgical dressings
- Skin graft sites
- Treatment of first degree burns
- Abrasions
- Wound packing

Item#	Description	Packaging
Sterile		
212	1" x 8"	50/bx
213	5" x 9"	50/bx
214	2" x 2"	25/bx
215	4" x 4"	25/bx
216	4" x 3 yds	6/bx

Connect with us!

Find our full line of products online at:
www.dukal.com

Telephone: 1-631-656-3800
Toll Free: 1-800-243-0741
Fax: 1-631-656-3810
Email: customerservice@dukal.com
To Place Orders: orders@dukal.com

Corporate Office

DUKAL Corporation
2 Fleetwood Court,
Ronkonkoma, NY 11779

