



PASSION. PARTNERSHIP. POSSIBILITIES.

TECHMED[®] Adjustable Aluminum Cane

Description

Easy push button adjustment for height settings



Features & Benefits

- Sure lock secures handle at any height adjustment
- Height adjusts in 1" increments (29"-38")
- Regular cane has latex free hand grip
- Quad cane has latex free hand grip made of soft hypalon and chrome plated steel base

Item#	Description	Packaging
6830-6	Regular Base	6/cs
6820-4	5" x 7" Small Quad Base	4/cs

Connect with us!

Find our full line of products online at:

www.dukal.com

Telephone: 1-631-656-3800

Toll Free: 1-800-243-0741

Fax: 1-631-656-3810

Email: customerservice@dukal.com

To Place Orders: orders@dukal.com

Corporate Office

DUKAL Corporation

2 Fleetwood Court,

Ronkonkoma, NY 11779

