



PASSION. PARTNERSHIP. POSSIBILITIES.

TECHMED® Foot Stools

Description

11" x 14" platform topped with 10" x 13" non-skid rubber mat



Features & Benefits

- Standard: Available in a chrome-plated or black finish, 450 lb. weight capacity
- Bariatric: Heavy duty chrome plated 16-gauge steel, 700 lb. weight capacity
- Handrail height 34"
- 9" from ground to platform
- Platform is 11" x 14"
- Not made with natural rubber latex

Item#	Description	Packaging
4351	Standard, Chrome	1/ea
4351BLK	Standard, Black	1/ea
4349	Standard, Handrail, Chrome	1/ea
4350	Bariatric, Chrome	1/ea
4348	Bariatric, Handrail, Chrome	1/ea

Connect with us!

Find our full line of products online at:

www.dukal.com

Telephone: 1-631-656-3800

Toll Free: 1-800-243-0741

Fax: 1-631-656-3810

Email: customerservice@dukal.com

To Place Orders: orders@dukal.com

Corporate Office

DUKAL Corporation

2 Fleetwood Court,

Ronkonkoma, NY 11779



Control # SL070855