



PASSION. PARTNERSHIP. POSSIBILITIES.

# TECHMED® IV Stands

## Description

Highly polished chrome plated steel



## Features & Benefits

- Adjust from 47" to 84"
- Locking collar allows pole to be set at desired height
- Four caster base has diameter of 23" and 15 lb weight capacity
- Five caster base has diameter of 20" and 22 lb weight capacity

| Item# | Description         | Packaging |
|-------|---------------------|-----------|
| 4353  | Four Caster, 2 Hook | 1/ea      |
| 4355  | Five Caster, 2 Hook | 1/ea      |
| 4356  | Five Caster, 4 Hook | 1/ea      |

### Connect with us!

Find our full line of products online at:

[www.dukal.com](http://www.dukal.com)

Telephone: 1-631-656-3800

Toll Free: 1-800-243-0741

Fax: 1-631-656-3810

Email: [customerservice@dukal.com](mailto:customerservice@dukal.com)

To Place Orders: [orders@dukal.com](mailto:orders@dukal.com)

### Corporate Office

DUKAL Corporation

2 Fleetwood Court,

Ronkonkoma, NY 11779

