



PASSION. PARTNERSHIP. POSSIBILITIES.



# Suture Removal Kit

## Description

Single-use kit that includes everything needed for safe and sanitary removal of sutures

## Features & Benefits

- Equipped with
  - 1 - Littauer scissors
  - 1 - Plastic forceps
  - 1 - Gauze pad, 3" x 3"
- Sterile with an easy to open peel back lid
- Not made with natural rubber latex

Item#	Description	Packaging
15100	Suture Removal Kit	50 trays/cs

### Connect with us!

Find our full line of products online at: [www.dukal.com](http://www.dukal.com)

Telephone: 1-631-656-3800  
Toll Free: 1-800-243-0741  
Fax: 1-631-656-3810  
Email: [customerservice@dukal.com](mailto:customerservice@dukal.com)  
To Place Orders: [orders@dukal.com](mailto:orders@dukal.com)

### Corporate Office

DUKAL Corporation  
2 Fleetwood Court,  
Ronkonkoma, NY 11779

