



PASSION. PARTNERSHIP. POSSIBILITIES.



# Barrier Film

## Description

DUKAL Barrier Film is a universal barrier with a low tack adhesive backing that is perfect for use in areas that are difficult to clean, disinfect or sterilize. Perforated 4" x 6" sheets are easy to apply, easy to remove and leave no residue. "Non-stick" edge facilitates removal. Packaged in a dispenser box.

## Features & Benefits

- Light tack adhesive allows for easy placement and removal with no residue left behind
- Perforated 4" x 6" sheets for desired fit
- Easy to use and disposable
- Covers most surface areas
- Available in clear or blue film

Item#	Description	Packaging
<b>Non-Sterile</b>		
27604	Barrier film - Clear	1200 sheets /rl, 12 rl/cs
27605	Barrier film - Blue	1200 sheets /rl, 12 rl/cs

### Connect with us!

Find our full line of products online at: [www.dukal.com](http://www.dukal.com)

Telephone: 1-631-656-3800  
 Toll Free: 1-800-243-0741  
 Fax: 1-631-656-3810  
 Email: [customerservice@dukal.com](mailto:customerservice@dukal.com)  
 To Place Orders: [orders@dukal.com](mailto:orders@dukal.com)

### Corporate Office

DUKAL Corporation  
 2 Fleetwood Court,  
 Ronkonkoma, NY 11779

