



PASSION. PARTNERSHIP. POSSIBILITIES.



# Sterilization Pouches

## Description

DUKAL Self-Sealing Sterilization Pouches are constructed with dual indicators (internal and external) that change color when completely processed by ETO or steam sterilization. These indicators provide a visual conformation that the optimum sterilization conditions were reached inside the sterilization unit and inside the pouch. Fabricated of rugged medical-grade paper and plastic film, these pouches guard against tears and perforations while providing the ease of use that every office is looking for.

## Features & Benefits

- Self-sealing pouch has a perforated tab for easy folding and to help maintain a perfect seal
- Triple sealed with extra wide “chevron” to protect instrument from piercing pouch
- Thumb notch for easy opening
- Blue see-through film allows for easy visibility inside pouch
- Rugged medical-grade paper guards against tears and perforations
- For use with steam and ETO sterilization only
- Lead free and not made with natural rubber latex

Item#	Description	Packaging
<b>Non-Sterile</b>		
27301	Sterilization Pouches, 2.25" x 4"	200/bx, 20 bx/cs
27302	Sterilization Pouches, 2.75" x 10"	200/bx, 20 bx/cs
27303	Sterilization Pouches, 3.5" x 5.25"	200/bx, 30 bx/cs
27304	Sterilization Pouches, 3.5" x 10"	200/bx, 20 bx/cs
27305	Sterilization Pouches, 3.5" x 10"	500/bx, 4 bx/cs
27306	Sterilization Pouches, 5.25" x 10"	200/bx, 15 bx/cs
27307	Sterilization Pouches, 7.5" x 13"	200/bx, 5 bx/cs

### Connect with us!

Find our full line of products online at: [www.dukal.com](http://www.dukal.com)

Telephone: 1-631-656-3800  
Toll Free: 1-800-243-0741  
Fax: 1-631-656-3810  
Email: [customerservice@dukal.com](mailto:customerservice@dukal.com)  
To Place Orders: [orders@dukal.com](mailto:orders@dukal.com)

### Corporate Office

DUKAL Corporation  
2 Fleetwood Court,  
Ronkonkoma, NY 11779

