



PASSION. PARTNERSHIP. POSSIBILITIES.

Level 2 Isolation Gowns

Features & Benefits

- Made of multi-ply SMS material
- Proprietary over-the-head design with easy tear away removal
- Elastic wrist with thumb hooks to provide a more secure fit under gloves
- Meets AAMI Level 2 requirements for fluid protection
- Used in low risk situations

Item#	Description	Packaging
Non-Sterile		
309	Universal - Yellow	10/bg; 10 bg/case
309XL	X-large - Yellow	10/bg; 10 bg/case



Connect with us!

Find our full line of products online at: www.dukal.com

Telephone: 1-631-656-3800
Toll Free: 1-800-243-0741
Fax: 1-631-656-3810
Email: customerservice@dukal.com
To Place Orders: orders@dukal.com

Corporate Office

DUKAL Corporation
2 Fleetwood Court,
Ronkonkoma, NY 11779



Effective Date: 7/31/20
ART-00053-DOC REV 1