



PASSION. PARTNERSHIP. POSSIBILITIES.



Digital Thermometer

Description

Quick 10 Second read with LCD display and low battery indicator.

Same Product,
New Quantity!
Now sold by the case

Features & Benefits

- Reads both Centigrade and Fahrenheit; Range: 32° – 43.9°C, 90° – 109.9° F
- Alarm/Beeper when reading is complete
- Battery life: 200 hours of continuous operation
- Auto shut off after 9 minutes of non use
- Built in memory enabling display of the previously measured value
- Water resistant

| Item# | Description | Packaging | Unit of Measure |
|-----------|----------------------------|--------------------|-----------------|
| 3210D-360 | Dual Scale, 10 second read | 12 ea/bx, 30 bx/cs | Case |

Connect with us!

Find our full line of products online at: www.dukal.com

Telephone: 1-631-656-3800
Toll Free: 1-800-243-0741
Fax: 1-631-656-3810
Email: customerservice@dukal.com
To Place Orders: orders@dukal.com

Corporate Office

DUKAL Corporation
2 Fleetwood Court,
Ronkonkoma, NY 11779



Effective Date: 07/13/20
ART-00049-DOC